



ASELight 50

LED EXAMINATION LIGHT

Unique coverage of varied and highest lighting requirements now available in one single compact casing.

asetronics 
Advanced Swiss Electronics



Application

- » General medical examinations and treatment in medical practices
- » Dermatological examinations
- » Wound care and minor surgical procedures
- » Dental treatment

Advantages of the new examination light in LED technology

- » Light intensity: 50,000 lux @ 100 cm distance
- » Colour temperature: 4'000 K (daylight quality)
- » Colour rendering index (CRI): 97
- » Light field diameter: 28 cm
- » Changeable light field for dental applications
- » Can be dimmed on the panel or for easy handling the integrated contact-free sensor can be used
- » No infrared radiation
- » High reliability and no maintenance costs
- » LED service life of over 70,000 hours
- » Extremely low energy consumption (typ. < 23 W) and extremely low energy costs (30 % in comparison to halogen)
- » Environmentally friendly (recyclable, does not contain any heavy metals)
- » High shock and vibration resistance
- » No risk of implosion (compared to halogen)
- » Safe because of the low voltage
- » Simple and quick to clean and disinfect thanks to high-quality materials and a smooth, sealed surface
- » Up-to-date design and technology
- » 5 years warranty on the light head

Swiss made

Installation:



Mobile stand



Mobile stand flex



Table



Ceiling mount



Wall mount

asetronics 

Advanced Swiss Electronics

Freiburgstrasse 251
3018 Bern
Switzerland

Internet: www.aselight.ch
www.asetronics.ch
Telephone: +41 (0)31 329 32 19
Telefax: +41 (0)31 329 31 22
E-mail: info@asetronics.ch

CE



(Version 09/2021)